

## Model DB-03

# Adjustable Doctor Blade (可調式刮刀)

Adjustable film applicators are designed to be convenient and easy to clean.

(可調式刮刀設計成方便且容易清潔)

The DB-03 instruments gap can be adjust from 10um to 3500um in 10um increments .

(Standard coating width : 100mm)



## SPECIFICATION(規格)

- Coating Width(塗佈寬度):
  - 50mm
  - 100mm
  - 150mm
  - 200mm
  - 250mm
  - 300mm
- Material(材質): Stainless (不鏽鋼)
- Gap (間隙範圍) : 10um to 3500um
- Gap interval (間隙間隔): 10um

## Operation notice (操作说明)

1. 歸零：刮刀置於平坦的(大理石)平台，轉動上端兩個分厘卡，使刮刀底部與兩邊的側板底部在同一平面，此時兩個分厘卡的刻度應該相同。分厘卡上的“0” 應對下端的準刻度線。
2. 同步旋轉上端的兩個分厘卡，調整到所需要的塗佈間隙，也就是塗佈厚度
3. 放置刮刀至塗佈機的塗佈平台
4. 倒入漿料後啟動塗佈機進行塗佈
5. 塗佈後以酒精或其他適當溶劑清潔刮刀，並保持於乾燥環境避免生鏽
6. 使用過程中需小心取放刮刀，否則經碰撞後將無法保證精度或使用